



ARES MEDİKAL SANAYİ TİCARET LTD. ŞTİ.



PRODUCT CATALOGUE

ARES

Ares Medikal Sanayi Ticaret Limited Şirketi founded in 2004, located in İzmir city of Turkey.



Ares Medikal wants to create easy access, economic, high quality products in market. we have a limited energy for produce your specific and new product needsings.

Ares Medikal is the leader company in his region about exporting medical products to many countries.


Ares Medikal has CE and ISO13485 certificates and follows all international quality rules.



CE 0068



Introducer Needle



Ares Introducer Needles are available in variable configurations and sizes ensuring that you have a needle specific to your procedural requirements.

Ergonomic hub ensures the perfect handling while the stainless steel needle shaft is particular for safe puncture.



INTRODUCER NEEDLE

CATALOG NO	Product	Size	Length
A2950	Introducer Needle Straight Type	18G	7 cm
A1090	Introducer Needle w/Wings	18G	7 cm
A3612	Transradial Introducer Needle	21G	4 cm



Radial

Femoral

Introducer Sets



Our expertise in precision brings the dilator and sheath transition to an excellent level.

Mindfully designed valve assures reliable haemostasis and facilitates catheter manipulation remarkably.

The smooth transition between components and dual tapered design minimize the risk of trauma in surrounding tissue of puncture point.

FEMORAL INTRODUCER SETS

CATALOG NO	Product	Size	Length
A3087	Introducer Set	4F	11 cm
A1250	Introducer Set	5F	11 cm
A1229	Introducer Set	6F	11 cm
A1236	Introducer Set	7F	11 cm
A1243	Introducer Set	8F	11 cm
A2677	Introducer Set	9F	11 cm
A2684	Introducer Set	10F	11 cm
A3384	Introducer Set	6F	23 cm
A3391	Introducer Set	7F	23 cm
A3407	Introducer Set	8F	23 cm

RADIAL INTRODUCER SETS

CATALOG NO	Product	Size	Length
A4992	Introducer Set	4F	7 cm
A2608	Introducer Set	5F	7 cm
A2615	Introducer Set	6F	7 cm
A5005	Introducer Set	5F	11 cm
A5012	Introducer Set	6F	11 cm





Radial



Femoral

PTFE Guidewires



CATALOG NO	Product	Diameter	Length	Tip Shape
A2837	Guide Wire	0,032"	150 cm	J-Type
A2868	Guide Wire	0,032"	180 cm	J-Type
A2844	Guide Wire	0,032"	260 cm	J-Type
A1120	Guide Wire	0,035"	150 cm	J-Type
A1151	Guide Wire	0,035"	180 cm	J-Type
A1137	Guide Wire	0,035"	260 cm	J-Type
A1168	Guide Wire	0,038"	150 cm	J-Type
A1199	Guide Wire	0,038"	180 cm	J-Type
A1175	Guide Wire	0,038"	260 cm	J-Type
A2851	Guide Wire	0,032"	150 cm	Straight
A1144	Guide Wire	0,035"	150 cm	Straight
A1182	Guide Wire	0,038"	150 cm	Straight

Ares guidewires are made of excellent quality torque-stable cores with optimal maneuverability and enhanced kink resistance particular for effective and atraumatic navigation through challenging vascular anatomies.

Guidewires are available with straight, angled and j-curved tip configurations. The proximal shaft with the ideal stiffness assures perfect pushability while the distal tip enables the operator for smooth and precise advancement.




Hydrophilic Guidewires



Innovative coating - Advanced coating technology produces a layer of moisture between the wire surface and either vessel walls or instrumentation

- Promotes smooth, frictionless navigation and maintains optimal tactile feel
- Maintains lubricity through multiple passes with minimal swelling

CATALOG NO	Product	Diameter	Length	Tip Shape
A2875	Guide Wire	0,035"	150 cm	Angled
A2899	Guide Wire	0,035"	180 cm	Angled
A2905	Guide Wire	0,035"	260 cm	Angled
A2912	Guide Wire	0,035"	150 cm	Straight




Manifold

Increased usability with smoothly rotating handles is particular to meet all procedural needs



High Pressure Stopcock



CATALOG NO	Product
A1069	Manifold Off Right
A3841	Manifold Off Left
A1076	Manifold Set
A1014	Angiography Syringe 12 ml

CATALOG NO	Product
A1113	500 PSI High Pressure Stopcock


Pressure: 500 PSI, holding higher pressure than most competitors
 Innovative ergonomic handles - smooth turning torque design improves ease of use and convenience.

Elevated handle design enables turning of the handles without touching the ports

High profile design for easy cleaning and disinfecting


Transparent housing for visual control

Available with 2 or 3 ports, OFF or ON handles, Right or Left orientation, different handle colors including color coded handles,




Hemostasis Valve Set

Easy manipulation
Absolute hemostasis
Maximum efficiency



Hemostasis Valve Set



Push-Pull Y Connector



Double Y Connector



Push-Click Y Connector



Screw Y Connector

Easy manipulation, absolute hemostasis with maximum efficiency

The highest level hemostatic valve with a large inner lumen to allow safe insertion of various interventional devices.

Enabling one-hand, one-touch control and absolute operational success.

Crystal clear body for perfect visualization


Side port with extension line for sound connection with other accessories.

CATALOG NO	Product
A1281	Screw Hemostasis Valve
A1274	Push-Pull Hemostasis Valve
A1380	Double Hemostasis Valve
A1298	Hemostasis Valve-Screw (Full Set)
A1304	Hemostasis Valve-Push-Pull (Full Set)
A5104	Hemostasis Valve Push-Click
A5128	Hemostasis Valve Push-Click (Full Set)

Inflation Device



Effortless Inflation
Successful Operation



 **Sidra**
One-Hand One-Touch



Classic
Perceivable Power



Gun Model
Remarkable Design

Effort-Free Error-Free


- One-hand use for easy and quick inflation and deflation
- The rotating handle and ergonomic device body is particular for comfortable grip and handling
- Features a 30 ATM pressure rating with a 30 mL syringe volume.
- The pressure gauge is available with 35 ATM option.



CATALOG NO	Product
A1052	Inflation Device 20ml, 30atm Classic Type
A2967	Inflation Device 30ml, 30atm Classic Type
A3780	Inflation Device 20ml, 30atm Gun Type
A3797	Inflation Device 30ml, 30atm Ergonomic Type
A2967	Ptca Kit (Balloon Inflation Device + Hemostasis Valve Set)


FEMORAL & RADIAL

Angio Closure Pad



Femoral Ped

Radial Ped Screw Type



Radial Ped w/Syringe

Luer lock fits all types of syringes

Air balloon with valve system provides wonderful pressure control

Large and middle sizes are available

CATALOG NO	Product
A1212	Radial Closure Pad w/Syringe
A2981	Femoral Closure Pad
A3827	Screw Type Radial Closure Pad



Radial



Femoral

Hemodialysis Catheter Kits




CONTENTS:

- 1 pcs .038"X70 cm. J / Flex nitinol guidewire
- 1 pcs Catheter
- 2 pcs Injection caps
- 2 pcs 10F+12F Dilators
- 1 pcs 18GA Needle
- 1 pcs Scalpel
- 1 pcs 5 ml Syringe
- 1 pcs Sticky Compress Pad

CATALOG NO	Lumen Type	Size	Model	Length
A3056	Double Lumen	11,5F	Curved	16 cm
A3049	Double Lumen	11,5F	Straight	16 cm
A3070	Double Lumen	11,5F	Curved	20 cm
A3063	Double Lumen	11,5F	Straight	20 cm
A3018	Double Lumen	12F	Curved	16 cm
A3001	Double Lumen	12F	Straight	16 cm
A3032	Double Lumen	12F	Curved	20 cm
A3025	Double Lumen	12F	Straight	20 cm



Connecting Tubing



PRESSURE LINE 800 PSI

CATALOG NO		Length
M/F	M/M	
A1519	A1571	30 cm
A1526	A1588	60 cm
A1533	A1595	90 cm
A1540	A1601	120 cm
A1557	A1618	150 cm
A1564	A1625	200 cm

PRESSURE LINE 1200 PSI

CATALOG NO		Length
FIXED MALE CONNECTOR	ROTATABLE MALE CONNECTOR	
A1632	A1694	30 cm
A1649	A1700	60 cm
A1656	A1717	90 cm
A2226	A2233	100 cm
A1663	A1724	120 cm
A1670	A1731	150 cm
A1687	A1748	200 cm


PRESSURE LINE 1200 PSI, ROTATING ADAPTER, BRAIDED TUBE

CATALOG NO		Length
VALVE ADAPTER	METALLIC ADAPTER	
A1755	A4329	30 cm
A1762	A4336	60 cm
A1779	A4343	90 cm
A2240	A4350	100 cm
A1789	A4367	120 cm
A1793	A4374	150 cm
A1809	A4381	200 cm




FEATURES

- Excellent flexibility
- Contains no plasticizers, phthalates or latex additives
- Custom profiles and tolerances
- Custom length and pressure rates
- Meets ISO 10993 guidelines for body contact applications
- CE approved
- Consistent performance
- Braided tube options
- Available in 84sh hardness



Connecting Tubing



PRESSURE LINE, SINGLE VALVE

CATALOG NO	Length
A2158	35 cm
A2165	120 cm
A2172	150 cm


PATIENT LINE, w/CHECK VALVES

CATALOG NO	Length
A2189	150 cm
A2202	250 cm
A3643	320 cm



PRESSURE LINE, COILED

CATALOG NO	Length
A1908	120 cm
A1915	150 cm
A1922	180 cm
A1939	200 cm
A1946	250 cm



CONNECTING TUBING, TWO WAY, w/CHECK VALVES

CATALOG NO		Length
Y - CONNECTOR	T - CONNECTOR	
A2097	A1977	150 cm
A2103	A1984	200 cm
A2110	A1991	250 cm
A2127	A2004	300 cm



GAS SAMPLING LINE

CATALOG NO	Product	Length
A4541	Gas Sampling Line	250 cm
A1960	Gas Sampling Line	300 cm

DSA MRI and CT Syringes






CATALOG NO	INJECTOR MODEL	INTERNAL VOLUME ML	ACCESSORY	DEVICE TYPE	UNIT
CT-200-M1	MEDTRON ACCUTRON CT	200	1PCS SYRINGE, EXTENSION LINE, J-TUBE	SINGLE HEAD	TOMOGRAPHY
CT-200-M1/Dual	MEDTRON ACCUTRON CT	200 / 200	2 PCS SYRINGE, 2-WAY EXTENSION LINE, SPIKE, J-TUBE	DOUBLE HEAD	TOMOGRAPHY
CT-190-A1	IMAXEON SALIENT	190	1PCS SYRINGE, EXTENSION LINE, J-TUBE	SINGLE HEAD	TOMOGRAPHY
CT-200-A1	MEDRAD VÍSTRON CT	200	1PCS SYRINGE, EXTENSION LINE, J-TUBE	SINGLE HEAD	TOMOGRAPHY
CT-200-A1	SINOPOWER S	200	1PCS SYRINGE, EXTENSION LINE, J-TUBE	SINGLE HEAD	TOMOGRAPHY
CT-200-A2	MEDRAD STELLANT	200	1PCS SYRINGE, EXTENSION LINE, J-TUBE	SINGLE HEAD	TOMOGRAPHY
CT-200-A2/Dual	MEDRAD STELLANT	200 / 200	2 PCS SYRINGE, 2-WAY EXTENSION LINE, SPIKE, J-TUBE	DOUBLE HEAD	TOMOGRAPHY
CT-200-B1	LF CT9000	200	1PCS SYRINGE, EXTENSION LINE, J-TUBE	SINGLE HEAD	TOMOGRAPHY
CT-200-B1/Y	LF OPTIVANTAGE	200	1 PCS SYRINGE, 2-WAY EXTENSION LINE, SPIKE, J-TUBE	DOUBLE HEAD	TOMOGRAPHY
CT-200-NE	NEMOTO	200	1PCS SYRINGE, EXTENSION LINE, J-TUBE	SINGLE HEAD	TOMOGRAPHY
CT-200-B1/Dual	LF OPTIVANTAGE	200 / 200	2 PCS SYRINGE, 2-WAY EXTENSION LINE, SPIKE, J-TUBE	DOUBLE HEAD	TOMOGRAPHY
CT/DSA-150-B1	LF ANGIOMAT ILLUMENA	150	1PCS SYRINGE, EXTENSION LINE, J-TUBE	SINGLE HEAD	TOMOGRAPHY / DSA
CT/DSA-200-B1	LF ANGIOMAT ILLUMENA	200	1PCS SYRINGE, EXTENSION LINE, J-TUBE	SINGLE HEAD	TOMOGRAPHY / DSA
DSA-150-B1	LF ANGIOMAT 6000	150	1PCS SYRINGE, J-TUBE	SINGLE HEAD	DSA
DSA-150-A1	MEDRAD MARK-V PROVIS	150	1PCS SYRINGE, J-TUBE	SINGLE HEAD	DSA
DSA-200-A1	MEDRAD MARK-V PROVIS	200	1PCS SYRINGE, J-TUBE	SINGLE HEAD	DSA
DSA-150-A2	MEDRAD MARK-7 ARTERION	150	1PCS SYRINGE, J-TUBE	SINGLE HEAD	DSA
MRI-115-A1	MEDRAD SPECTRIS / SOLARIS	115	1PCS SYRINGE, EXTENSION LINE, J-TUBE	SINGLE HEAD	MR
MRI-65-A1	MEDRAD SPECTRIS / SOLARIS	65	1PCS SYRINGE, EXTENSION LINE, J-TUBE	SINGLE HEAD	MR
MRI-65/115-A1/Dual	MEDRAD SPECTRIS / SOLARIS	65 / 115	2 PCS SYRINGE, 2-WAY EXTENSION LINE, SPIKE, J-TUBE	DOUBLE HEAD	MR
MRI-60/60-B1/Dual	LF OPTISTAR	60 / 60	2 PCS SYRINGE, 2-WAY EXTENSION LINE, SPIKE, J-TUBE	DOUBLE HEAD	MR

Angiography Drape Set



CONTENT	SIZE (cm)	PIECES
Angiography Drape	150 x 240 cm	1
Tools Table Drape	100 x 150 cm	1
Spanch	7,5 x 7,5 cm	10
Solution Cup	500 cc	1
Fluoroscopy Cover	85 x 85 cm	1
Standart Surgery Gown	Large	2
Towel	40 x 40 cm	2
Pack Of Sterilization	100 x 100 cm	1






Guiding Catheter

304 Stainless steel braiding gives excellent stability and ultra endurance

Large inner lumen design assures outstanding contrast flow and device compatibility


Maximum guidewire capacity 0.070" (1.8mm)

Maximum pressure rate is 1200 psi (84 bar)





Curve	Type	5Fr		6Fr		7Fr		
		Catalogue No.	I.D.	Catalogue No.	I.D.	Catalogue No.	I.D.	
RIGHT	JUDKINS	JR3.0	GCJR305F	0.0578"	GCJR306F	0.0708"	GCJR307F	0.0807"
		JR3.5	GCJR355F	0.0578"	GCJR356F	0.0708"	GCJR357F	0.0807"
		JR4.0	GCJR405F	0.0578"	GCJR406F	0.0708"	GCJR407F	0.0807"
		JR4.5	GCJR455F	0.0578"	GCJR456F	0.0708"	GCJR457F	0.0807"
		JR5.0	GCJR505F	0.0578"	GCJR506F	0.0708"	GCJR507F	0.0807"
	AMPLATZ	AR1.0	GCAR105F	0.0578"	GCAR106F	0.0708"	GCAR107F	0.0807"
		AR2.0	GCAR205F	0.0578"	GCAR206F	0.0708"	GCAR207F	0.0807"
		AL0.75	GCAL755F	0.0578"	GCAL756F	0.0708"	GCAL757F	0.0807"
		AL1.0	GCAL105F	0.0578"	GCAL106F	0.0708"	GCAL107F	0.0807"
		AL1.5	GCAL155F	0.0578"	GCAL156F	0.0708"	GCAL157F	0.0807"
LEFT	JUDKINS	AL2.0	GCAL205F	0.0578"	GCAL206F	0.0708"	GCAL207F	0.0807"
		AL2.5	GCAL255F	0.0578"	GCAL256F	0.0708"	GCAL257F	0.0807"
		AL3.0	GCAL305F	0.0578"	GCAL306F	0.0708"	GCAL307F	0.0807"
		JL3.0	GCJL305F	0.0578"	GCJL306F	0.0708"	GCJL307F	0.0807"
		JL3.5	GCJL355F	0.0578"	GCJL356F	0.0708"	GCJL357F	0.0807"
	EXTRA BACKUP	JL4.0	GCJL405F	0.0578"	GCJL406F	0.0708"	GCJL407F	0.0807"
		JL4.5	GCJL455F	0.0578"	GCJL456F	0.0708"	GCJL457F	0.0807"
		JL5.0	GCJL505F	0.0578"	GCJL506F	0.0708"	GCJL507F	0.0807"
		EBU3.0	GCEBU305F	0.0578"	GCEBU306F	0.0708"	GCEBU307F	0.0807"
		EBU3.5	GCEBU355F	0.0578"	GCEBU356F	0.0708"	GCEBU357F	0.0807"
	JUDKINS	EBU3.75	GCEBU3755F	0.0578"	GCEBU3756F	0.0708"	GCEBU3757F	0.0807"
		EBU4.0	GCEBU405F	0.0578"	GCEBU406F	0.0708"	GCEBU407F	0.0807"
		EBU4.5	GCEBU455F	0.0578"	GCEBU456F	0.0708"	GCEBU457F	0.0807"
		EBU5.0	GCEBU505F	0.0578"	GCEBU506F	0.0708"	GCEBU507F	0.0807"

Angiography Catheter



Soft tip reduces potential for vessel trauma
Double braided stainless steel for 1:1 Torque control
Maximum Guide wire: 0.038" (0.97 mm.)
Maximum Pressure: 1200psi (84 bar)



	TIP SHAPE	Spesifikasyon	CATALOG NO
ANGIOGRAPHY CATHETER	AC-JL	CATETER CARD JUDKINS PERFORMA 4FR X 100CM 038 JL3.5	ACJL354F
		CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JL3.5	ACJL355F
		CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JL3.5	ACJL356F
		CATETER CARD JUDKINS PERFORMA 4FR X 100CM 038 JL4.0	ACJL404F
		CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JL4.0	ACJL405F
		CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JL4.0	ACJL406F
		CATETER CARD JUDKINS HIGHFLOW 5FR X 100CM 038 JL4.0	HACJL405F
		CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JL4.5	ACJL455F
		CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JL4.5	ACJL456F
		CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JL5.0	ACJL505F
		CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JL5.0	ACJL506F
		CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JL6.0	ACJL605F
		AC-JR	CATETER CARD JUDKINS PERFORMA 4FR X 100CM 038 JR3.5
	CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JR3.5		ACJR355F
	CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JR3.5		ACJR356F
	CATETER CARD JUDKINS PERFORMA 4FR X 100CM 038 JR4.0		ACJR404F
	CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JR4.0		ACJR405F
	CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JR4.0		ACJR406F
	CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JR5.0		ACJR505F
	CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JR5.0		ACJR506F
	CATETER CARD JUDKINS PERFORMA 5FR X 100CM 038 JR6.0		ACJR605F
	CATETER CARD JUDKINS PERFORMA 6FR X 100CM 038 JR6.0		ACJR606F
	AC-AL	CATETER CARD AMPLATZ PERFORMA 5FR X 100CM 038 AL1	ACAL15F
		CATETER CARD AMPLATZ PERFORMA 6FR X 100CM 038 AL1	ACAL16F
		CATETER CARD AMPLATZ PERFORMA 5FR X 100CM 038 AL2	ACAL25F
		CATETER CARD AMPLATZ PERFORMA 6FR X 100CM 038 AL2	ACAL26F
		CATETER CARD AMPLATZ PERFORMA 5FR X 100CM 038 AL3	ACAL35F
		CATETER CARD AMPLATZ PERFORMA 6FR X 100CM 038 AL3	ACAL36F
		CATETER CARD AMPLATZ PERFORMA 5FR X 100CM 038 AL4	ACAL45F
		CATETER CARD AMPLATZ PERFORMA 6FR X 100CM 038 AL4	ACAL46F
	AC-AR	CATETER CARD AMPLATZ PERFORMA 5FR X 100CM 038 AR1	ACAR15F
		CATETER CARD AMPLATZ PERFORMA 6FR X 100CM 038 AR1	ACAR16F
		CATETER CARD AMPLATZ PERFORMA 5FR X 100CM 038 AR2	ACAR25F
		CATETER CARD AMPLATZ PERFORMA 6FR X 100CM 038 AR2	ACAR26F
		CATETER CARD AMPLATZ PERFORMA 5FR X 100CM 038 AR3	ACAR35F
		CATETER CARD AMPLATZ PERFORMA 6FR X 100CM 038 AR3	ACAR36F
	AC-TIG	CATETER CARDIOLOGICO PERF ULTIMATE 1 RADIAL 5FR 100CM 038	ACTIG15F
		CATETER CARDIOLOGICO PERF ULTIMATE 2 RADIAL 5FR 100CM 038	ACTIG25F
		CATETER CARDIOLOGICO PERF ULTIMATE 1 RADIAL 6FR 100CM 038	ACTIG16F
	AC-PIG	CATETER PIGTAIL ULTRA HIGH FLOW 5FR X 65CM	ACPIGS5F1065
		CATETER CARD PIGTAIL 145° 6FR X 110CM 038	ACPIGA6F
		CATETER CARD PIGTAIL 145° 5FR X 110CM 038	ACPIGASF
		CATETER CARD PIGTAIL RETO 4FR X 80CM 038	ACPIGS4F1080
		CATETER CARD PIGTAIL RETO 4FR X 110CM 038	ACPIGS4F
		CATETER CARD PIGTAIL RETO 5FR X 110CM 038	ACPIG5F)
		CATETER CARD PIGTAIL RETO 5FR X 100CM 038 20 MARKER BANDS	ACPIG5F20M10H
		CATETER CARD PIGTAIL RETO 5FR X 125CM 038	ACPIGS5F
		CATETER CARD PIGTAIL RETO 6FR X 110CM 038	ACPIG6F)
		CATETER PIT TALL 5FR X 100CM 038 20BANDS(1CM APART) 5 SH	ACPIGT2015F
		CATETER PIT TALL 5FR X 100CM 038 20BANDS(1CM APART) 10 SH	ACPIGT2025F
		CATETER PIT TALL 5FR X 65CM 038 11BANDS(1CM APART) 10 SH	ACPIGT1115F165
		CATETER PIT TALL 5FR X 100CM 038 11BANDS(1CM APART) 10 SH	ACPIGT1115F
	CATETER PIT TALL 5FR X 65CM 038 2 MARKER BANDS(2CM APART) 10 SH	ACPIGT215F165	
	AC-LCB	CATETER CORONARIO BYPASS LEFT 5FR X 100CM 038 LCB	ACLCB5F
		CATETER CORONARIO BYPASS RIGHT 6FR X 100CM 038 RCB	ACRCB6F
		CATETER CORONARIO BYPASS RIGHT 5FR X 100CM 038 RCB	ACRCB5F
	AC-IM	CATETER MAMARIA INTERNA PERFORMA 5FR X 100 CM	ACIM5F
	AC-CIM	CATETER MAMARIA INTERNA PERFORMA 4FR X 100 CM 038	ACIM4F
		CATETER MAMARIA INTERNA PERFORMA 6FR X 100 CM	ACIM6F
	AC-MPA	CATETER CARD MULTIPURPOSE 4FR X 100CM 038 A2	ACMPA24F
		CATETER CARD MULTIPURPOSE 5FR X 100CM 038 B2	ACMPA25F
		CATETER CARD MULTIPURPOSE 4FR X 100CM 038 A1	ACMPA14F
		CATETER CARD MULTIPURPOSE 6FR X 100CM 038 A1	ACMPA16F
		CATETER CARD MULTIPURPOSE 5FR X 100CM 038 A1	ACMPA15F
		CATETER CARD MULTIPURPOSE 6FR X 100CM 038 A2	ACMPA26F



ARES MEDİKAL SANAYİ TİCARET LTD. ŞTİ.

Aşık Veysel Mah. 5821/1 Sk. No:6
Karabağlar 35110
İzmir, Turkey

Phone : +90 232 264 70 00
Fax : +90 232 264 90 00
Email : info@ares-medikal.com
Web : www.ares-medikal.com

