



GEOTEK[®]

Healthcare Products

ABOUT US / HAKKIMIZDA

Geotek is a manufacturer company of interventional medical devices, dedicated to continuous improvement of the performance and quality of its service and products.

Our aim is expanding our distribution around the world with new and improved products and we commit to take care of the customer's needs with manufacturing high quality medical device.

Milestones of the company

2003 was founded in Ankara, Turkey. Started with representations of well-known international medical companies and maintained a sub-dealers sales and marketing network all around Turkey in a short period

2009 opened production facility and got CE certification

2012 launched first Automatic Biopsy Needle

2017 got FDA 510k approval for Disposable Needle Guides

2020 New Reusable Gun Series Launched

2021 USA Branch Office Opened

2022 New Factory in Operation

Geotek sürekli kalite ve performans gelişimi ile dünya sağlık sektörüne tıbbi malzeme üretmektedir. Amacımız, dünya genelinde dağıtım ağıımızı, en yüksek kalitede geliştirilmiş ürünler ile geliştirerek, kullanıcıların beklentilerini karşılamaktır.

Geotek'in Kilometre Taşları;

2003 yılında kurulmuş, önceleri çok iyi bilinen markaların ithalatı ile başlamış, kısa sürede yurt içinde geniş bayii ağına ulaşmıştır

2009 yılında üretim bazlı fabrika kurmuş ve CE kalite belgesi almıştır

2012 yılında kendi adıyla ilk Otomatik Biyopsi İğnesi satışa sunulmuştur

2017 yılında Tek Kullanımlık İğne Kılavuzu için FDA onayı almıştır

2020 yılında Yeni Reusable Biyopsi Tabanca Serileri satışa sunulmuştur

2021 yılında ABD ofisi açılmıştır

2022 Yeni fabrika hizmete girmiştir



GEOTEK®
Healthcare Products

CE 1984

CONTENTS / İÇİNDEKİLER

UROLOGY CATHETERS / ÜROLOJİ KATATERLERİ

Long Durable DJ Ureteral Stent Set / <i>Karbotan DJ Üreteral Stent Set</i>	4
Double J Ureteral Stent Set / <i>Double J Üreteral Stent Set</i>	5
Ureteral Catheter / <i>Üreter Kateter</i>	6
Ureteral Balloon Dilator / <i>Üreteral Balon Dilatör</i>	7
Re-entry Malecot Catheter / <i>Re-Entry Malekot Kateter</i>	8
Amplatz Renal Dilator Set / <i>Amplatz Renal Dilatör Set</i>	9
Direct Puncture Nephrostomy Set / <i>Direkt Giriş Nefrostomi Seti</i>	10
Single J Stent Set / <i>Single J Üreteral Stent Set</i>	11
Suprapubic Drainage Catheter Set / <i>Suprapubik Drenaj Seti</i>	12
Percutaneous Nephrostomy Set / <i>Perkütan Nefrostomi Seti</i>	13
Stone Basket / <i>Taş Basket</i>	14
Stone Grasper / <i>Taş Tutucu</i>	15
Nephrostomy Balloon Dilator / <i>Nefrostomi Balon Kateter</i>	16
Ureteral Access Sheath / <i>Üreteral Giriş Kılıfı</i>	17

UROLOGY CATHETERS


ÜROLOJİ KATATERLERİ

LONG DURABLE DJ URETERAL STENT SET

KARBOTAN DJ ÜRETERAL STENT SET

- ↳ Made by Carbothane
- ↳ Suitable for long term implantation
- ↳ Higher kink resistance
- ↳ High radiopacity
- ↳ Depth marker
- ↳ Long pusher option between 40-100 cm
- ↳ Suture on distal tip

Diameter (F)	Length (cm)	Ref. No: (One End Open)	Ref. No: (Both End Open)
3	12	GDJ312-OEC-CAR	GDJ312-BEO-CAR
	14	GDJ314-OEC-CAR	GDJ314-BEO-CAR
	16	GDJ316-OEC-CAR	GDJ316-BEO-CAR
	20	GDJ320-OEC-CAR	GDJ320-BEO-CAR
4	12	GDJ412-OEC-CAR	GDJ412-BEO-CAR
	14	GDJ414-OEC-CAR	GDJ414-BEO-CAR
	16	GDJ416-OEC-CAR	GDJ416-BEO-CAR
	20	GDJ420-OEC-CAR	GDJ420-BEO-CAR
4.7	22	GDJ422-OEC-CAR	GDJ422-BEO-CAR
	24	GDJ424-OEC-CAR	GDJ424-BEO-CAR
	26	GDJ426-OEC-CAR	GDJ426-BEO-CAR
	12	GDJ4712-OEC-CAR	GDJ4712-BEO-CAR
4.8	16	GDJ4716-OEC-CAR	GDJ4716-BEO-CAR
	20	GDJ4720-OEC-CAR	GDJ4720-BEO-CAR
	22	GDJ4722-OEC-CAR	GDJ4722-BEO-CAR
	24	GDJ4724-OEC-CAR	GDJ4724-BEO-CAR
5	26	GDJ4726-OEC-CAR	GDJ4726-BEO-CAR
	28	GDJ4728-OEC-CAR	GDJ4728-BEO-CAR
	12	GDJ4812-OEC-CAR	GDJ4812-BEO-CAR
	20	GDJ4820-OEC-CAR	GDJ4820-BEO-CAR
6	22	GDJ4822-OEC-CAR	GDJ4822-BEO-CAR
	24	GDJ4824-OEC-CAR	GDJ4824-BEO-CAR
	26	GDJ4826-OEC-CAR	GDJ4826-BEO-CAR
	28	GDJ4828-OEC-CAR	GDJ4828-BEO-CAR
7	16	GDJ516-OEC-CAR	GDJ516-BEO-CAR
	22	GDJ522-OEC-CAR	GDJ522-BEO-CAR
	24	GDJ524-OEC-CAR	GDJ524-BEO-CAR
	26	GDJ526-OEC-CAR	GDJ526-BEO-CAR
8	28	GDJ528-OEC-CAR	GDJ528-BEO-CAR
	20	GDJ620-OEC-CAR	GDJ620-BEO-CAR
	22	GDJ622-OEC-CAR	GDJ622-BEO-CAR
	24	GDJ624-OEC-CAR	GDJ624-BEO-CAR
8	26	GDJ626-OEC-CAR	GDJ626-BEO-CAR
	28	GDJ628-OEC-CAR	GDJ628-BEO-CAR
	30	GDJ630-OEC-CAR	GDJ630-BEO-CAR
	24	GDJ724-OEC-CAR	GDJ724-BEO-CAR
8	26	GDJ726-OEC-CAR	GDJ726-BEO-CAR
	28	GDJ728-OEC-CAR	GDJ728-BEO-CAR
	30	GDJ730-OEC-CAR	GDJ730-BEO-CAR
	22	GDJ822-OEC-CAR	GDJ822-BEO-CAR
8	24	GDJ824-OEC-CAR	GDJ824-BEO-CAR
	26	GDJ826-OEC-CAR	GDJ826-BEO-CAR
	28	GDJ828-OEC-CAR	GDJ828-BEO-CAR
	30	GDJ830-OEC-CAR	GDJ830-BEO-CAR



- ↳ Karbotan materyal
- ↳ Uzun dönem kalıcı stent
- ↳ Yüksek bükülme direnci
- ↳ Yüksek radyopasite
- ↳ Derinlik ölççeği
- ↳ 40-100 cm arasında uzun itici seçeneği
- ↳ Radyoopak itici seçeneği
- ↳ Alt uçta sütün

↳ Set includes DJ stent + Pusher + Guide wire + Two Clamps

DOUBLE J URETERAL STENT SET


DOUBLE J ÜRETERAL STENT SET

- ↳ Soft polyurethane material
- ↳ Higher kink resistance
- ↳ High radiopacity
- ↳ Depth marker
- ↳ Long pusher option between 40 cm - 100 cm
- ↳ Suture on distal tip



Diameter (F)	Length (cm)	Ref. No: (One End Open)	Ref. No: (Both End Open)	
3	12	GDJ312-OEC	GDJ312-BEO	
	14	GDJ314-OEC	GDJ314-BEO	
	16	GDJ316-OEC	GDJ316-BEO	
	20	GDJ320-OEC	GDJ320-BEO	
4	12	GDJ412-OEC	GDJ412-BEO	
	14	GDJ414-OEC	GDJ414-BEO	
	16	GDJ416-OEC	GDJ416-BEO	
	20	GDJ420-OEC	GDJ420-BEO	
	22	GDJ422-OEC	GDJ422-BEO	
	24	GDJ424-OEC	GDJ424-BEO	
4.7	26	GDJ426-OEC	GDJ426-BEO	
	12	GDJ4712-OEC	GDJ4712-BEO	
	16	GDJ4716-OEC	GDJ4716-BEO	
	20	GDJ4720-OEC	GDJ4720-BEO	
	22	GDJ4722-OEC	GDJ4722-BEO	
	24	GDJ4724-OEC	GDJ4724-BEO	
Adjustable	26	GDJ4726-OEC	GDJ4726-BEO	
	28	GDJ4728-OEC	GDJ4728-BEO	
	22/30	VDJ47-OEC	VDJ47-BEO	
	4.8	12	GDJ4812-OEC	GDJ4812-BEO
		16	GDJ4816-OEC	GDJ4816-BEO
		20	GDJ4820-OEC	GDJ4820-BEO
22		GDJ4822-OEC	GDJ4822-BEO	
24		GDJ4824-OEC	GDJ4824-BEO	
26		GDJ4826-OEC	GDJ4826-BEO	
Adjustable	28	GDJ4828-OEC	GDJ4828-BEO	
	22/30	VDJ48-OEC	VDJ48-BEO	
	5	16	GDJ516-OEC	GDJ516-BEO
		20	GDJ520-OEC	GDJ520-BEO
		22	GDJ522-OEC	GDJ522-BEO
		24	GDJ524-OEC	GDJ524-BEO
26		GDJ526-OEC	GDJ526-BEO	
28		GDJ528-OEC	GDJ528-BEO	
6	20	GDJ620-OEC	GDJ620-BEO	
	22	GDJ622-OEC	GDJ622-BEO	
	24	GDJ624-OEC	GDJ624-BEO	
	26	GDJ626-OEC	GDJ626-BEO	
	28	GDJ628-OEC	GDJ628-BEO	
	30	GDJ630-OEC	GDJ630-BEO	
Adjustable	22/30	VDJ06-OEC	VDJ06-BEO	
	7	22	GDJ722-OEC	GDJ722-BEO
		24	GDJ724-OEC	GDJ724-BEO
		26	GDJ726-OEC	GDJ726-BEO
		28	GDJ728-OEC	GDJ728-BEO
		30	GDJ730-OEC	GDJ730-BEO
Adjustable		22/30	VDJ07-OEC	VDJ07-BEO
	8	22	GDJ822-OEC	GDJ822-BEO
		24	GDJ824-OEC	GDJ824-BEO
		26	GDJ826-OEC	GDJ826-BEO
		28	GDJ828-OEC	GDJ828-BEO
		30	GDJ830-OEC	GDJ830-BEO

- ↳ Yumuşak poliüretan materyal
- ↳ Yüksek bükülme direnci
- ↳ Yüksek radyopasite
- ↳ Derinlik ölçeği
- ↳ 40-100 cm arasında uzun itici seçeneği
- ↳ Radyopak itici seçeneği
- ↳ Alt uçta suture



DJ + PUSHER ONLY

Please Use Ref. No as GDJPXX-OEC or BEO


HYDROPHILIC DJ STENT SET

Please Use Ref. No as GDJXX-HOEC or HBEO

- ↳ Set includes DJ stent + Pusher + Guide wire + Two Clamps

URETERAL CATHETER

ÜRETERAL KATETER



- ↳ Wide size range
- ↳ Depth marker
- ↳ Luer lock adapter
- ↳ Metal mandren

- ↳ Geniş ölçü alternatifi
- ↳ Derinlik işareti
- ↳ Kilitli adaptör
- ↳ Metal mandren

URETERAL CATHETER CLOSED TIP

Ref. No	Diameter (F)	Length (cm)
GUC370-CE	3	70
GUC470-CE	4	70
GUC570-CE	5	70
GUC670-CE	6	70
GUC770-CE	7	70

URETERAL CATHETER OPEN TIP

Ref. No	Diameter (F)	Length (cm)
GUC370-OE	3	70
GUC470-OE	4	70
GUC570-OE	5	70
GUC670-OE	6	70
GUC770-OE	7	70

URETERAL BALLOON DILATOR

ÜRETERAL BALON DİLATÖR



- ↳ Used for dilatation of ureter
- ↳ 2 ways catheter
- ↳ Platinum markers at both end
- ↳ 18 to 20 atm rated burst pressure
- ↳ Supplied with stopcock

- ↳ Üreterin dilatasyonu için kullanılır
- ↳ 2 yönlü kateter shaftı
- ↳ Her iki ucunda platinyum işaret
- ↳ 18-20 atm max. basınç
- ↳ 3 yollu musluk

Ref. No	Balloon Diameter (mm)	Balloon Length (cm)	Shaft (F)	Usable Length (cm)	Recommended Guide wire (inch)	Balloon RBP (atm)
GBD0544	4	4	5	70	.38"	20
GBD0546	4	6	5	70	.38"	20
GBD5410	4	10	5	70	.38"	20
GBD0554	5	4	5	70	.38"	20
GBD0556	5	6	5	70	.38"	20
GBD5510	5	10	5	70	.38"	18
GBD0564	6	4	5	70	.38"	18
GBD0566	6	6	5	70	.38"	18
GBD5610	6	10	5	70	.38"	18
GBD0444	4	4	4	70	.18"	20
GBD0454	5	4	4	70	.18"	20

RE-ENTRY MALECOT CATHETER

RE-ENTRY MALEKOT KATETER

- ↳ Two flexible stylets
- ↳ High radiopacity
- ↳ Smooth surface
- ↳ Drainage bag adaptor

- ↳ Esnek stile
- ↳ Yüksek radyoopasite
- ↳ Girişimi kolaylaştıran yüzey
- ↳ Bağlantı adaptörü

Ref. No	Diameter (F)	Body Length (cm)	Tail Length (cm)
GRM 1430	14	30	15
GRM 1630	16	30	15
GRM 1830	18	30	15
GRM 2030	20	30	15

AMPLATZ RENAL DILATOR SET

AMPLATZ RENAL DİLATÖR SET

- ↳ Used for renal dilatation before stone removal procedures
- ↳ High radiopacity
- ↳ Wide range of dilator options: 6F to 30F
- ↳ 4 Pcs sheath and 70 cm introducer catheter




- ↳ Perkütan taş çıkarma işlemi öncesi renal dilatasyon için uygundur
- ↳ Yüksek radyoopasite
- ↳ Geniş dilatör seçeneği 6F-30F
- ↳ 4 Adet kılıf ve 70 cm introducer kateter

Ref. No	Description	Size
GUDS-12	Amplatz Renal Dilator Set	6 F - 30 F
GUDS-22	Amplatz Renal Dilator w/PTFE Sheath	22 F
GUDS-24	Amplatz Renal Dilator w/PTFE Sheath	24 F
GUDS-26	Amplatz Renal Dilator w/PTFE Sheath	26 F
GUDS-28	Amplatz Renal Dilator w/PTFE Sheath	28 F
GUDS-30	Amplatz Renal Dilator w/PTFE Sheath	30 F

DIRECT PUNCTURE NEPHROSTOMY SET


DİREKT GİRİŞ NEFROSTOMİ SETİ

- ↳ Pigtail nephrostomy catheter with metal stiffener and trocar point stylet
- ↳ Easy penetration with sharp beveled trocar and tapered pigtail catheter tip
- ↳ Large drainage holes
- ↳ Hydrophilic coated and locked catheter types are optional



Ref. No	Description	Size (Fr)	Length (cm)
GDPN08	Direct Puncture Nephrostomy Set	8	30
GDPN10	Direct Puncture Nephrostomy Set	10	30
GDPN12	Direct Puncture Nephrostomy Set	12	30
GDPN14	Direct Puncture Nephrostomy Set	14	30
GDPN08-L	Direct Puncture Nephrostomy Set with Locking	8	30
GDPN10-L	Direct Puncture Nephrostomy Set with Locking	10	30
GDPN12-L	Direct Puncture Nephrostomy Set with Locking	12	30
GDPN14-L	Direct Puncture Nephrostomy Set with Locking	14	30
GDPNH08	Direct Puncture Nephrostomy Set Hydrophilic Coated	8	30
GDPNH10	Direct Puncture Nephrostomy Set Hydrophilic Coated	10	30
GDPNH12	Direct Puncture Nephrostomy Set Hydrophilic Coated	12	30
GDPNH14	Direct Puncture Nephrostomy Set Hydrophilic Coated	14	30
GDPNH08-L	Direct Puncture Nephrostomy Set with Locking and Hydrophilic Coted	8	30
GDPNH10-L	Direct Puncture Nephrostomy Set with Locking and Hydrophilic Coted	10	30
GDPNH12-L	Direct Puncture Nephrostomy Set with Locking and Hydrophilic Coted	12	30
GDPNH14-L	Direct Puncture Nephrostomy Set with Locking and Hydrophilic Coted	14	30

- ↳ Pigtail nefrostomi kateteri ve içine yerleştirilmiş trokar uçlu giriş iğnesinden oluşur
- ↳ Keskin uçlu giriş iğnesi ve pigtail kateterin özgün uç yapısı sayesinde kolay girişim
- ↳ Büyük drenaj delikleri
- ↳ Hidrofilik kaplı kateter ve kilit sistemi seçenekleri



SINGLE J URETERAL STENT SET

SINGLE J ÜRETER STENT SET

- ↳ Two different colors Single J catheters for left and right side
- ↳ Guide wire
- ↳ Dilators 7 F to 9 F
- ↳ Urine Bag Connector

- ↳ Single J Kateter
- ↳ Kılavuz tel
- ↳ Dilatör 7 FR - 9 FR
- ↳ İdrar Torba Bağlantısı




Ref. No	Diameter (F)	Length (cm)
GSJ 670	6	70
GSJ 770	7	70
GSJ 870	8	70

SUPRAPUBIC DRAINAGE CATHETER SET

SUPRAPUBİK DRENAJ SETİ

- ↳ %100 Silicone Foley Catheter
- ↳ Trocar with peel away sheath
- ↳ Scalpel




- ↳ %100 Silikon Foley Kateter
- ↳ Yırılabilir kılıflı Trokar
- ↳ Bisturi

Ref. No	Catheter Diameter (F)	Catheter Length (cm)
GSP10	10	35
GSP12	12	35
GSP14	14	35
GSP16	16	35
GSP18	18	35

PERCUTANEOUS NEPHROSTOMY SET

PERKÜTAN NEFROSTOMİ SETİ

- ↳ Pigtail or malecot tip catheter with:
- ↳ Puncture needle
- ↳ Dilators
- ↳ Stiff guide wire or lunderquist guide wire options
- ↳ Scalpel
- ↳ Drainage bag connector



- ↳ Pigtail Nefrostomi kateter ve:
- ↳ Girişim iğnesi
- ↳ Dilatörler
- ↳ Kılavuz tel
- ↳ Bistüri
- ↳ Drenaj torba bağlantı adaptörü

NEPHROSTOMY SET

Ref. No (Malecot)	Ref. No (Pigtail)	Diameter (F)	With Locking	Hydrophilic	Hydrophilic & Locking
-	GPNS06	6	GPNS06-L	GPNSH06	GPNSH06-L
-	GPNS08	8	GPNS08-L	GPNSH08	GPNSH08-L
GMNS10	GPNS10	10	GPNS10-L	GPNSH10	GPNSH10-L
GMNS12	GPNS12	12	GPNS12-L	GPNSH12	GPNSH12-L
GMNS14	GPNS14	14	GPNS14-L	GPNSH14	GPNSH14-L
GMNS16	GPNS16	16	GPNS16-L	GPNSH16	GPNSH16-L

NEPHROSTOMY CATHETER

Ref. No (Malecot)	Ref. No (Pigtail)	Diameter (F)	With Locking	Hydrophilic	Hydrophilic & Locking
-	GNC06	6	GNC06-L	GNCH06	GNCH06-L
-	GNC08	8	GNC08-L	GNCH08	GNCH08-L
GMC10	GNC10	10	GNC10-L	GNCH10	GNCH10-L
GMC12	GNC12	12	GNC12-L	GNCH12	GNCH12-L
GMC14	GNC14	14	GNC14-L	GNCH14	GNCH14-L
GMC16	GNC16	16	GNC16-L	GNCH16	GNCH16-L

STONE BASKET

TAŞ BASKET

- ↳ Stainless steel and nitinol alternatives
- ↳ Designed for endoscopic manipulation
- ↳ Detachable handle

- ↳ Paslanmaz çelik ve nitinol seçeneği
- ↳ Endoskopik kullanıma uygun dizayn
- ↳ Ayrılabilir sap



STAINLESS STEEL

Ref. No	Diameter (F)	Length (cm)	BASKET TYPE	WIRE TYPE
GSBH390	3	90	HELICAL	Round
GSBS390			PARALLEL	Flat
GSBH3110	3	110	HELICAL	Round
GSBS3110			PARALLEL	Flat
GSBH490	4	90	HELICAL	Round
GSBS490			PARALLEL	Flat
GSBH4110	4	110	HELICAL	Round
GSBS4110			PARALLEL	Flat
GSBH590	5	90	HELICAL	Round
GSBS590			PARALLEL	Flat
GSBH5110	5	110	HELICAL	Round
GSBS5110			PARALLEL	Flat

NITINOL

Ref. No	Diameter (F)	Length (cm)	BASKET TYPE	WIRE TYPE	Zero Tip Ref. No
GSBN1990	1.9	90	HELICAL	Round	GSBZ1990
GSBN19110	1.9	110	HELICAL	Round	GSBZ19110
GSBN19120	1.9	120	HELICAL	Round	GSBZ19120
GSBN2290	2.2	90	HELICAL	Round	GSBZ2290
GSBN22110	2.2	110	HELICAL	Round	GSBZ22110
GSBN390	3	90	HELICAL	Round	GSBZ390
GSBN3110	3	110	HELICAL	Round	GSBZ3110
GSBN490	4	90	HELICAL	Round	GSBZ490
GSBN4110	4	110	HELICAL	Round	GSBZ4110
GSBN590	5	90	HELICAL	Round	GSBZ590
GSBN5110	5	110	HELICAL	Round	GSBZ5110

STONE GRASPER


TAŞ TUTUCU

- ↳ Stainless steel
- ↳ Designed for endoscopic manipulation
- ↳ Detachable handle

- ↳ Paslanmaz çelik
- ↳ Endoskopik kullanıma uygun ayrılabilir dizayn
- ↳ Ayrılabilir sap



Ref. No	Description	Diameter (F)	Length (cm)
GSG390	Grasping Forceps, Three Pong	3	90
GSG490	Grasping Forceps, Three Pong	4	90
GSG590	Grasping Forceps, Three Pong	5	90



NEPHROSTOMY BALLOON DILATOR

NEFROSTOMİ BALON KATETER

- ↳ Used for dilatation of renal tract
- ↳ 2 ways catheter
- ↳ Platinum markers at both end
- ↳ 17 atm rated burst pressure
- ↳ Supplied with 30 F PTFE Sheath

- ↳ Perkütan taş prosedürlerinde renal dilatasyon için kullanılır
- ↳ İki yollu kateter vardır
- ↳ Balonun her iki ucunda Platinyum marker
- ↳ 17atm max. basınç
- ↳ 30F PTFE kılıfla birlikte sağlanır



Ref. No	Balloon Diameter (mm)	Balloon Length (cm)	Shaft (F)	Usable Length (cm)	Recommended Guide wire (inch)	Balloon RBP (atm)
GNBD10	10	12	7	55	.38"	17



URETERAL ACCESS SHEATH

ÜRETERAL GİRİŞ KILIFI




↳ Ureteral access sheath provides ureteral dilation and continuous working channel for the introduction of endoscopes and instruments during ureteral access procedures

↳ Üreteral giriş prosedürü süresince dilatasyona, endoskop ve enstrümanların girişi için sürekli bir çalışma kanalına imkan tanır

Ref. No	Diameter (F)	Length (cm)
GUA1125	9.5/11.5	25
GUA1138	9.5/11.5	38
GUA1145	9.5/11.5	45
GUA1225	10/12	25
GUA1238	10/12	38
GUA1245	10/12	45
GUA1425	12/14	25
GUA1438	12/14	38
GUA1445	12/14	45





NOTES



A series of horizontal lines for writing, spanning the width of the page. The lines are evenly spaced and extend across most of the page width, leaving a margin on the left. The background features a light gray grid pattern and a large, faint yellow circle that overlaps the writing area.

NOTES



GEOTEK®

Healthcare Products

Head Office: İvedik OSB Ağaç-Metal San. Sit.
1436. Sokak No: 12 06370 Ostim - Ankara / TÜRKİYE
Tel: +90 (312) 395 98 52 • Fax: +90 (312) 395 98 53

Factory: Eskişehir Yolu 42. Km ASO 2 ALCI OSB
2014. cadde No: 11 06909 Sincan - Ankara / TÜRKİYE
Tel: +90 (312) 514 64 04

e-mail: info@geotekmedical.com

www.geotekmedical.com

USA Office: 333 5th Avenue 5th Floor, New York, NY 10016, USA

Tel: +1 (908) 699-0662 • info@geotekmedicalusa.com

www.geotekmedicalusa.com

